






# Standardized voltage settings

For better resolution and experimental reproducibility

# Why standardized settings are needed?

## Lack of reproducibility in science



Baker, M. 1,500 scientists lift the lid on reproducibility. *Nature* **533**, 452–454 (2016).

# What is responsible for irreproducible findings in cytometry experiments?


<b>Statistics</b>	Sample size, significance cutoffs
<b>Reagents</b>	Lot-to-lot variability, titration, specificity, reagent interactions
<b>Instruments</b>	Detector settings, optical filters, compensation, fluidic instability
<b>Humans</b>	Operator-to-operator inconsistency
<b>Analysis</b>	Data cleanup, manual gating vs computational approach

# Different approaches to setting voltages

	By eye	QC	Copy from template	Matching target MFIs
Fast?	-	+++	+++	++
Reproducible?	-	+++	+	+++
Best separation?	?	+	?	+++


# How standardized settings are determined

For each fluorescence detector: record cells stained with an appropriate anti-CD4 conjugate over range of PMT voltages



# How standardized settings are determined

Calculate stain index and plot SI vs voltage. Determine EC90 voltage.



## Regression Model Options

Regression Type

FourParmLogistic


Name	Fixed	Last Used
Imin	0	0
Imax		52.754
IC50		435.09
HillSlope		9.4433
r2		0.996476798829

Number of regression dots

30



# Standardized settings don't depend on the cell type or the actual fluorochrome used



Jones, DD, et. al., U Penn Flow Cytometry & Cell Sorting Resource Laboratory



# MFI values should stay consistent, not voltages

- Beads were recorded with the EC90 voltages to obtain “target” MFI values
- Before each experiment, the same beads will be run and voltage adjusted to obtain the same MFI value



# How to use standardized settings

## 1. Create new experiment using the “RFCC new experiment” template

The screenshot displays a flow cytometry software interface with a Global Worksheet titled "MFI target matching" and an Inspector window for "Cytometer Settings".

**Global Worksheet - MFI target matching**

The worksheet contains a grid of plots and summary tables. The plots show MFI target matching for various parameters, including beads, singlets, BV421 peak 7, BV510 peak 7, BV605 peak 8, PE peak 7, BV650 peak 8, PE-Cy5 peak 7, BV711 peak 8, PE-Cy5-5 peak 7, BV786 peak 8, FITC peak 8, APC peak 6, PerCP peak 7, Alexa 700 peak 6, BV650 peak 8, PE-Cy5 peak 7, BV711 peak 8, PE-Cy5-5 peak 7, BV396 378/29-A, and APC-Cy7 peak 8.

**Summary Tables:**

Parameter	Population	Median
FITC-A	23,035	####
BV421-A	35,879	####
PE-A	38,295	####
APC-A	29,120	####
PerCP-Cy5-5-A	23,285	####
BV510-A	21,308	####
PE-CF594-A	51,050	####
APC-Alexa 700-A	33,137	####
BUV396 378/29-A	1,151	####
BV605-A	49,745	####
PE-Cy5-A	52,406	####
APC-Cy7-A	21,455	####
BUV737 740/35 6...	10,345	####
BV650-A	29,919	####
PE-Cy5-5-A	25,811	####
BV711-A	23,414	####
PE-Cy7-A	29,792	####
BV786-A	12,686	####

**Inspector - Cytometer Settings**

The Inspector window shows the Cytometer Settings for the experiment. The Parameters tab is active, displaying a list of parameters and their corresponding settings.

Parameter	Voltage	...	A	H	W
FSC	350	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
SSC	214	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
FITC	480	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
PerCP-Cy5-5	572	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
BV421	336	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
BV510	362	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
BV605	436	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
BV650	439	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
BV711	481	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
BV786	488	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
PE	427	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
PE-CF594	502	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
PE-Cy5	435	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
PE-Cy5-5	445	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
PE-Cy7	510	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
APC	551	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
APC-Alexa 700	561	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
APC-Cy7	482	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
BUV396 378/29	370	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
BUV737 740/35 690LP	464	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

# How to use standardized settings

## 2. Delete unwanted parameters

The screenshot displays a flow cytometry software interface with a grid of histograms and a list of parameters. The histograms show 'MFI target matching-8-peak beads' for various parameters. The list of parameters includes:


Population	Parameter	Median	Value
<input checked="" type="checkbox"/>	FITC peak 8	FITC-A	23,035
<input checked="" type="checkbox"/>	BV421 peak 7	BV421-A	35,879
<input checked="" type="checkbox"/>	PE peak 7	PE-A	38,295
<input checked="" type="checkbox"/>	APC peak 6	APC-A	29,120
<input checked="" type="checkbox"/>	PerCP peak 7	PerCP-Cy5-5-A	23,285
<input checked="" type="checkbox"/>	BV510 peak 7	BV510-A	21,308
<input checked="" type="checkbox"/>	PE-CF594 peak 7	PE-CF594-A	51,050
<input checked="" type="checkbox"/>	Alexa 700 peak 6	APC-Alexa 700-A	33,137
<input checked="" type="checkbox"/>	BUV395 peak 6	BUV396 378/29-A	1,151
<input checked="" type="checkbox"/>	BV605 peak 8	BV605-A	49,745
<input checked="" type="checkbox"/>	PE-Cy5 peak 7	PE-Cy5-A	52,406
<input checked="" type="checkbox"/>	APC-Cy7 peak 8	APC-Cy7-A	21,455
<input checked="" type="checkbox"/>	BUV737 peak 8	BUV737 740/35 6...	10,345
<input checked="" type="checkbox"/>	BV650 peak 8	BV650-A	29,919
<input checked="" type="checkbox"/>	PE-Cy5-5 peak 7	PE-Cy5-5-A	25,811
<input checked="" type="checkbox"/>	BV711 peak 8	BV711-A	23,414
<input checked="" type="checkbox"/>	PE-Cy7 peak 8	PE-Cy7-A	29,792
<input checked="" type="checkbox"/>	BV786 peak 8	BV786-A	12,686

The 'Inspector - Cytometer Settings' window is open, showing a table of parameters and their settings:

Parameter	Voltage	Log	A	H	W
FSC	350	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
SSC	214	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
FITC	480	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
BV421	336	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
BV711	481	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
PE	427	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
APC-Cy7	408	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
APC	551	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

# How to use standardized settings

## 3. Run beads and adjust voltage to match target MFI






# How to use standardized settings

4. Create compensation controls and look at full stain:
  - Adjust FSC and SSC voltage for your cells
  - Confirm all populations are on scale
  - **Only** adjust voltage to get positive events on scale!
5. Record controls and calculate compensation as usual
6. Load analysis template to apply plots and gates from previous experiment

# Benefits to using standardized settings

- Best resolution
- Experimental rigor and reproducibility
  - Easier to compare data over time with fewer batch effects
  - Confidence that any change you see is biological
  - What if a reviewer asks how you set the voltages?
- Can quantify how Fortessa 2 is different from Fortessa 3
- Instrument specific information to aid in panel design
  - Fluorochrome brightness
  - Spread matrix

# Instrument specific fluorochrome brightness



Use for critical tertiary antigens  
where resolution is most important

Use for primary antigens to decrease  
spread: CD45, CD4, CD8, etc







# Instrument specific spread matrix

...will increase spread of this parameter

Staining with this dye...

	FITC 530/30	PerCP-Cy5-5 710/50	BV421 450/50	BV510 525/50	BV605 610/20	BV650 670/30	BV711 710/50	BV786 780/60	PE 586/15	PE-CF594 610/20	PE-Cy5 670/30	PE-Cy5-5 710/50	PE-Cy7 780/60	APC 670/30	APC-Alexa 700 710/50	APC-Cy7 780/60	BUV396 378/29	BUV737 740/35
FITC		34.8	16.3	17.4	14.7	14.1	16.8	17.9	11.9	11.9	26.6	25	13.6	18.5	35.9	13.6	21.2	11.4
PerCP-Cy5.5	14.1		16.8	16.3	14.7	27.7	45.6	41.3	11.4	11.9	47.8	48.9	31.5	53.8	113	25	21.7	23.4
BV421	15.2	29.9		20.7	16.8	14.7	16.8	17.9	11.9	13	27.2	25	14.1	19	36.9	14.1	22.3	11.4
BV510	15.8	34.2	17.9		46.2	39.7	35.3	32	20.1	22.8	27.2	25	13.6	28.2	58.2	14.1	21.2	17.4
BV605	14.7	72.3	18.5	17.4		77.2	71.2	56.5	66.3	76.7	88.1	52.2	31.5	52.7	106	14.7	21.2	28.2
BV650	15.2	68	19.6	17.4	48.4		94.7	70.7	22.3	29.9	84.8	50	28.8	108	180	26.6	22.3	33.7
BV711	15.2	144	21.7	16.3	15.7	33.7		137	11.9	11.9	45.1	61.4	49.4	62.5	365	64.7	22.3	63
BV786	14.7	35.9	25.5	17.4	15.2	16.8	30.4		11.9	11.9	27.7	25.5	33.1	22.8	62	38	21.7	28.2
PE	16.8	107	16.3	16.9	36.4	28.8	29.3	21.2		109	99.6	54.9	28.8	32	48.9	14.1	21.7	14.7
PE-CF594	14.7	201	16.8	16.3	45.6	39.7	45.1	27.2	79.4		177	99.6	59.8	47.8	69.6	15.2	21.7	18.5
PE-Cy5	15.2	694	20.1	16.9	23.4	76.7	131	51.1	33.1	26.6		237	137	267	321	52.7	22.8	42.4
NY700	15.2	286	16.8	16.9	15.7	33.1	60.3	23.9	50.5	37.5	212		89.2	212	802	86	21.7	22.8
PE-Cy7	15.2	58.7	17.9	15.8	14.1	14.7	20.1	82.7	28.8	20.6	33.7	35.8		20.6	44	38.6	21.7	22.8
APC	15.2	84.3	15.8	16.3	15.2	30.4	33.7	27.2	11.9	15.2	146	62	40.7		244	41.8	21.7	25
Alexa 700	15.2	42.9	16.8	16.3	13.6	13.6	27.2	28.2	11.4	11.9	32.6	41.3	33.7	31.5		41.8	21.2	22.3
APC-Cy7	14.7	35.9	16.8	16.3	13.6	14.7	18.5	45.6	11.4	11.9	42.9	29.3	90.3	53.2	77.2		21.2	18.5
BUV395	15.2	29.9	16.3	16.9	13.6	13	16.3	17.4	11.4	11.9	26.6	23.9	14.1	17.9	36.4	13.6		11.4
BUV737	15.2	166	17.9	16.9	14.1	22.8	49.4	60.9	11.9	11.9	31.5	52.2	50	33.1	238	68.5	22.8	



# Future goals

- Measure CD4 single stains and generate target MFIs for the following RFCC instruments
  - Fusion (done) and Ariall
  - FACSymphony S6 (coming soon)
  - Bigfoot
  - MA900
  - All 3 Fortessas
- Publish instrument specific spread matrix and fluorochrome brightness ranking on RFCC website